

OxyMask, OxyKid & OxyTyke Test

Please complete the test and forward to your Supervisor.

* Your Name:

* Hospital, Clinic, Company:

Department:

* Employee ID:

* Supervisor's e-mail:

* **1. What FiO₂ range can OxyMask and OxyKid deliver?**

24% - 50%

24% - 90%

40% - 70%

60% - 80%

2. At what flow rate can OxyMask and OxyKid achieve it's highest FiO₂?

10 LPM

15 LPM

20 LPM

40 LPM

3. What FiO₂ range can OxyTyke deliver?

22% - 50%

24% - 90%

24% - 50%

60% - 80%

4. What is the minimum flow rate achievable by the OxyTyke?

10 ML/M

15 ML/M

250 ML/M

40 ML/M

5. Name a clinical advantage of the OxyMask and OxyKid.

One mask can deliver the full range of O₂ Therapy

The open mask is less claustrophobic & well tolerated

FiO₂ Titration is simple by adjusting the flow meter up or down

All the above

6. How do you titrate O₂ concentration on the OxyMask, OxyKid and OxyTyke?

Replace the colored Venturi Jets to the specific FiO₂ required

Adjust Flow Rate up or down to maintain required SpO₂

You can't. You must change to another mask.