REF

OT-1025-8 - OxyMask™ Tyke
OT-7025-8 - OxyMask™ Tyke Non PVC

Intended Use: A minimal contact, open oxygen delivery device for use with USP medical grade oxygen.

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Caution: Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury or death. Do not smoke while using this device. Single use product not to be reused due to effect on performance, cross contamination and risk of infection.

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EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

CH REP

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

UK CA 0086 CE 2797

Emergo Consulting (UK) Limited
c/o Cr 360—UL International
Compass House, Vision Park Histon
Cambridge, CB24 9BZ
United Kingdom

Southmedic

50 Alliance Blvd.
Barrie, Ontario
Canada L4M 5K3
705-726-9383
705-728-9537 FAX
thebetteroxygenmask.com

EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

CH REP

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

UK CA 0086 CE 2797

Emergo Consulting (UK) Limited
c/o Cr 360—UL International
Compass House, Vision Park Histon
Cambridge, CB24 9BZ
United Kingdom

Oxygen Flow Settings	Approximate O ₂ Concentrations
0.25 LPM	22%
0.50 LPM	25%
0.75 LPM	30%
1 LPM	35%
1.5 LPM	43%
2 LPM	47%
3 LPM	50%
≥5 LPM	65%

*Actual FIO₂ will vary. Breath rate, depth of breathing and distance of the diffuser from the face will affect performance.

** Actual flow rate required should be prescribed in conjunction with pulse oximetry.

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