
OxyMask™

- OM-1125-8 - OxyMask™ Adult 7' Tubing
- OM-7025-8 - OxyMask™ Adult Non PVC 7' Tubing
- OM-1425-8 - OxyMask™ Adult Ear Elastics 7' Tubing
- OM-1125-14 - OxyMask™ Adult 14' Tubing
- OP-1125-8 - OxyMask™ Plus 7' Tubing
- OP-7025-8 - OxyMask™ Plus Non PVC 7' Tubing
- OP-1425-8 - OxyMask™ Plus Ear Elastics 7' Tubing
- OP-1125-14 - OxyMask™ Plus 14' Tubing
- OK-1125-8 - OxyMask™ Kid 7' Tubing
- OK-7025-8 - OxyMask™ Kid Non PVC 7' Tubing
- OK-1125-14 - OxyMask™ Kid 14' Tubing

REF

Intended Use: A minimal contact, open oxygen delivery device for use with USP medical grade oxygen.

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Caution: Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury or death. Do not smoke while using this device. Single use product not to be reused due to effect on performance, cross contamination and risk of infection.

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Made in Canada

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

IFU0288G

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Emergo Consulting (UK) Limited
c/o Cr 360—UL International
Compass House, Vision Park Histon
Cambridge, CB24 9BZ
United Kingdom

Made in Canada

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

IFU0288G

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Emergo Consulting (UK) Limited
c/o Cr 360—UL International
Compass House, Vision Park Histon
Cambridge, CB24 9BZ
United Kingdom

Oxy₂Mask™

thebetteroxygenmask.com

*Actual FIO₂ will vary. Breath rate, depth of breathing and distance of the diffuser from the face will affect performance.

** Actual flow rate required should be prescribed in conjunction with pulse oximetry.

Oxygen Flow Settings	Approximate O ₂ Concentrations
1 LPM	24% - 27%
2 LPM	27% - 32%
3 LPM	30% - 60%
4 LPM	33% - 65%
5 LPM	36% - 69%
7 LPM	48% - 80%
10 LPM	53% - 85%
12 LPM	57% - 89%
≥15 LPM	60% - 90%

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