Intended Use: A minimal contact, open oxygen delivery device for use with USP medical grade oxygen for monitoring breathing by providing a means to sample exhaled CO₂.

Warning: Hypercapnia will not be reliably detected by sampling through this device, the displayed CO₂ values or any other parameter/index driven from CO₂ will vary with the sampling port position, oxygen flow rate and respiratory pattern.

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Caution: Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury or death. Do not smoke while using this device. Single use product not to be reused due to effect on performance, cross contamination and risk of infection. The device has a 5-year shelf life from the date of manufacture. Date of manufacture is stated on labelling. Store device in original packaging until use. Inspect product and packaging prior to use. Do not use if product appears to be damaged.

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O₂ / EtCO₂

thebetteroxygenmask.com

Oxygen Flow Settings	Approximate O ₂ Concentrations
1 LPM	24% - 34%
2 LPM	30% - 50%
4 LPM	37% - 55%
6 LPM	47% - 70%
8 LPM	50% - 79%
10 LPM	54% - 80%
12-15 LPM	59% - 83%

*Actual FIO₂ will vary. Breath rate, depth of breathing and distance of the diffuser from the face will affect performance.
** Actual flow rate required should be prescribed in conjunction with pulse oximetry.

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REPORTING SUSPECTED SIDE EFFECTS

If you suspect that use of this medical device is associated with an adverse occurrence please report the details of the incident to your health authority, and your local distributor or directly to the manufacturer at www.southmedic.com.

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