

Finally ... A HIGHLY EFFICIENT AEROSOL MASK



OxyMask™

AEROSOL MASK
An OxyArm™ Product



OxyTyke Aerosol
OTN-1025



OxyKid Aerosol
OKN-1025



OxyMask Aerosol
OMN-1025



OxyPlus Aerosol
OPN-1025
15% larger than OxyMask Aerosol

OxyMask Aerosol Series



OxyMask Aerosol Series Logical Advantages:

- **Greater patient compliance**
because mask is more comfortable to wear
- **Swivel adapter**
 - allows for optimal aerosol delivery when patient is laying down or is vertical
 - tapered design allows interfacing with a wide variety of inner and outer elbow nebulizer cups
- **Unique design provides directional flow** to mouth and nose generating a higher percentage of respirable dose and greater aerosol efficiency
- **Closed sides and open bottom design:**
 - greater aerosol efficiency
 - prevents aerosol from going into the eyes
 - minimizes environmental contamination
 - minimizes exposure to caregiver and other patients



OxyTyke Aerosol KIT
OTN-4025-8
OTN-4025-8USC



OxyKid Aerosol KIT
OKN-4025-8
OKN-4025-8USC



OxyMask Aerosol KIT
OMN-4025-8
OMN-4025-8USC



OxyPlus Aerosol KIT
OPN-4025-8
OPN-4025-8USC
15% larger than
OxyMask Aerosol

Aerosol flow WITHOUT obstruction:



Competitor masks express aerosol as shown upward into the eye area

Note: Tested at 8 LPM




The OxyMask Aerosol Series express aerosol towards patient's nose and mouth

Aerosol flow WITH obstruction:



Competitor masks express aerosol as shown upward into the eye area



The OxyMask Aerosol Series express aerosol downward

Southmedic

705-726-9383
email: smedic@southmedic.com

www.southmedic.com www.oxyarm.com

50 Alliance Blvd., Barrie, ON., Canada L4M 5K3
ISO 13485

OxyMask™
AEROSOL MASK
An OxyArm™ Product



Emergo Europe
Molenstraat 15
2513 BH
The Hague
The Netherlands

Phone: +31(0) 70 345 8570
Fax: +31(0) 70 346 7299