


NOW ONE

Endotracheal Tube Protector - Eco Friendly WITH SECURITY RING

From I.D. 6.5 - 9.5 ADULT
From I.D. 3.5 - 6.0 PEDIATRIC



Note positioning for optimal
compression resistance and
placement of adjustable Tie-wrap

BITE BLOCK advantages:

- Bite Block's unique design fits multiple sized endotracheal tubes (ET Tubes)
- Designed to overlap on smaller ET Tubes, giving the device even greater compression resistance
- Soft enough to protect teeth and yet extremely durable to protect ET Tube
- Has Slip-Security Ring that fits over the standard ET Tube 15 mm adaptor
- Easy to attach and adjust
- Clear to visualize tube placement
- Adjustable Tie-wrap to secure Bite Block
- Material is fungus resistant
- Environmentally Friendly – Made from Non PVC material

ETPSR-00A50	(box of 50 Adult)
ETPSR-00A100	(box of 100 Adult)
ETPSR-00P50	(box of 50 Pediatric)
ETPSR-00P100	(box of 100 Pediatric)

Southmedic 

1-705-726-9383 1-705-728-9537 Fax
email: smedic@southmedic.com
www.southmedic.com

50 Alliance Blvd., Barrie, ON., Canada L4M 5K3
ISO 13485

P0029 07/09

EcoSmedic TM

BITE BLOCK
Made from Eco Friendly Non-PVC Material