

SPSS Metal Handle Safety Scalpel System

Safety Tab for Optimum Protection

Exclusive design effectively guards against accidental injury during installation.

X-Ray Detectable

Cartridge and Tab are X-Ray Detectable

Metal Handle

Handle has the same weight and feel of conventional scalpels

Reusable Handle

SPSS are more cost effective because they create less disposable waste.

Polymer Coated Surgical Blade

Regarded by many as the sharpest, smoothest cutting edge available.

Easy Loading and Disposal



The cartridge system provides a safe and easy way to use and dispose of blades.

Clear Plastic Sheath

Design incorporates thumb rest for easy manipulation.

Customized Ruler

Designed into reverse side of handle.



	Sterile 50/Box 3 Boxes/Case	Non Sterile 100/Bag 5 Bags/Case
10	73-8010	73-8210
11	73-8011	73-8211
15	73-8015	73-8215
15c	73-8015c	73-8215c
20	73-8020	73-8220
21	73-8021	73-8221
22	73-8022	73-8222
23	73-8023	73-8223
24	73-8024	73-8224

SPSS Safety Scalpel Metal Handles

73-8030-1	#3	1/pkg
73-8031-1	#3L	1/pkg
73-8040-1	#4	1/pkg
73-8041-1	#4L	1/pkg
73-8071-1	#7	1/pkg

THE BLADE YOU NEVER TOUCH

The Weight You Expect

- A Metal Handle Safety Scalpel System with the weight and balance of traditional scalpels
- SPSS cartridges are paired with Southmedic's Polymer Coated Blade that reduces incision drag and tissue trauma
- It's innovative features meet OSHA Needle Stick Safety & Prevention Guidelines
- Truly a blade you never touch

Southmedic 
www.southmedic.com **Blades • Scalpels**
Polymer Coated

SPSS Disposable Safety Scalpel System

The Blade you Never Touch - with Tear Away Safety Tab

- Fully retracts to expose entire blade
- SPSS cartridges are paired with Southmedic's Polymer Coated Blade that reduces tissue drag and trauma
- Its innovative features meet OSHA Needle Stick Safety & Prevention Guidelines
- The perfect "Disposable" Safety Scalpel System



	Sterile 10/Box 10 Boxes/Case	Non Sterile 100/Bag 5 Bags/Case
--	---	--

10	73-9010	73-9210
11	73-9011	73-9211
15	73-9015	73-9215
15c	73-9015c	73-9215c
20	73-9020	73-9220
21	73-9021	73-9221
22	73-9022	73-9222
23	73-9023	73-9223
24	73-9024	73-9224

Southmedic 
Blades + Scalpels
 Polymer Coated

1-800-463-7146 Toll Free in North America
 1-705-726-9383 1-705-728-9537 Fax
 www.southmedic.com email: smedic@southmedic.com
 50 Alliance Blvd., Barrie, ON., Canada L4M 5K3
 ISO 13485
 Made in Canada

