

DELIVERS FREEDOM

OxyChin™

An OxyArm™ Product
1st Open Oxygen System™
We are the experts!
Patent Pending

Directional
Diffuser
Technology

Safe

- Provides quick access for oral suctioning
- Ideal for conscious sedation & post op oxygenation
- When emesis occurs

Comfortable

- Open design allows communication between patient and caregiver
- Reduces claustrophobia
- Minimal contact oxygen delivery system

Efficient

- OxyChin delivers 1-15 LPM with one device
- Rate Trend Monitoring model delivers O₂ and measures CO₂

Oxygen Flow Settings	Approximate O ₂ Concentrations
1 LPM	21% - 27%
2 LPM	28% - 31%
4 LPM	32% - 35%
6 LPM	36% - 39%
8 LPM	40% - 43%
10 LPM	42% - 47%
12-15 LPM	48% - 50%

ORDER CODES:

O₂ MODELS

- OC-1025 - box of 25
- OC-1025-8 - with tubing - box of 25

CO₂ MODELS

- OC-2025 - Rate Trend Monitoring - box of 25
- OC-2025-8 - with tubing - Rate Trend Monitoring - box of 25

Southmedic

50 Alliance Blvd., Barrie, ON., Canada L4M 5K3
705-726-9383 • 705-726-9537 Fax ISO 13485 EDI 406000117
Email: smedic@southmedic.on.ca

www.oxyarm.com • www.southmedic.com

P0017A - May06