

Finally ... a simple an **MULTIPURPOSE MASK**

High Oxygen 24%-80% • High Humidity • Aerosol Therapy

OxyMulti Mask™

Oxygen, Humidity and Aerosol Therapy
An OxyArm™ Product

*Ideal for Acute Care / EMS
situations when you need
High FiO₂, High Humidity
and/or delivery of
Medicated Aerosol
Therapy*

Clinical Advantages:

- One Mask System
- FiO₂ = 24% to 80%
- Flow Rate 1-15 LPM
- High Humidity Mask
- Efficient Aerosol Delivery Mask
- Swivel Adaptor
 - OD is 22mm • ID is 15mm tapered



Finally ... a simple MULTIPURPOSE MASK

High Oxygen 24% - 80% • High Humidity • Aerosol Therapy

OxyMulti Mask

OHH-1025-8

OPHH-1025-8

OHH-1425-8

OPHH-1425-8

OMMK-4025-8

OPMK-4025-8

OMMK-1425-8

OPMK-1425-8

2174-PT-04

OxyMulti Mask – with Adaptor and 7' tubing

OxyMulti Plus Mask – with Adaptor and 7' tubing

OxyMulti Mask – with Adaptor, 7' Tubing and Ear Elastics


OxyMulti Plus Mask – with Adaptor, 7' Tubing and Ear Elastics

OxyMulti Mask Kit – with Adaptor, 7' Tubing, and Nebulizer

OxyMulti Plus Mask Kit – with Adaptor, 7' Tubing and Nebulizer

OxyMulti Mask Kit – with Adaptor, 7' Tubing, Nebulizer and Ear Elastics

OxyMulti Plus Mask Kit – with Adaptor, 7' Tubing, Nebulizer and Ear Elastics

15mm OD Tapered Adaptor 

OxyTrach

OAT-1025

OAT-4025-8

TCS-1005

OxyTrach Mask – High Oxygen/High Humidity/Aerosol Therapy/with Adaptor

OxyTrach Mask Kit – with Mask, 7' Tubing and Nebulizer/with Adaptor

Comfort Strap for OxyTrach Mask



OxyTrach

OxyTrach

Clinical Advantages:

- One Mask System
- FiO₂ = 24% to 80%
- Flow Rate 1-15 LPM
- High Humidity Mask
- Efficient Aerosol Delivery Mask
- Swivel Adaptor
 - OD is 22mm • ID is 15mm tapered

Ideal for Acute Care / EMS situations when you need High FiO₂, High Humidity and/or delivery of Medicated Aerosol Therapy



Comfort Strap TCS-1005

OxyMulti Mask

Oxygen, Humidity and Aerosol Therapy

An OxyArm™ Product

1st Open Oxygen System™

We are the experts!

Patent Pending

OXYMASK IS A REGISTERED TRADEMARK

Southmedic

705-726-9383

email: smedic@southmedic.com

www.southmedic.com www.oxyarm.com

50 Alliance Blvd., Barrie, ON., Canada L4M 5K3
ISO 13485

11/10 P0044