

# DELIVERING FREEDOM

OxyArm™



## Testimonials

"My quality of life has improved significantly. I wish I had it a year ago."  
George, Saskatoon, Canada



"I love my OxyArm. I take it to the pool to do my aqua fitness. My tank sits on the side of the pool and I am free to do my exercise."  
Boulder, Colorado, USA

"I feel accepted by my friends. They think it's cool and would like one too."  
Barrie, Ontario



"The light weight design is pleasant and comfortable. I have true confidence in this device as it has been the easiest to use for me to date. The light headband is also a pleasure"  
Jim, Barrie, Ontario

A PARADIGM SHIFT IN OXYGEN THERAPY

# OXYARM DELIVERS FREEDOM

IDEAL FOR BURN UNIT, TRAUMA CENTRE AND HOMECARE

## Diffuser Flow Rates

Oxygen Flow Settings	Approximate O <sub>2</sub> Concentrations
1 LPM	24% - 27%
2 LPM	28% - 31%
4 LPM	32% - 35%
6 LPM	36% - 39%
8 LPM	40% - 43%
10 LPM	42% - 47%
12-15 LPM	48% - 50%

## OxyCannula Flow Rates

Oxygen Flow Settings	1-8 LPM




Diffuser Arm O<sub>2</sub>




OxyCannula

- Adapts for left or right sided positioning
- Minimal contact design - no facial contact
- Designed for both nose & mouth breathers
- Relief from nasal irritation
- Appealing design - high tech, sleek and light weight
- Latex-free construction
- Delivers freedom





Product Codes  
**OA-PLUS-1125**  
**OA-PLUS-1125-8 (7ft tubing)**

## OxyArm™ - Adjustable version



Product Codes  
**OA-HB-1125**    **OA-DA-1025**  
**OA-HB-1100**    **OA-DA-1006**



Product Codes  
**OA-PLUS-3125**  
**OA-PLUS-3125-8**

## OxyArm™ - Adjustable version



Product Codes  
**OA-HB-1125**    **OA-DA-3025**  
**OA-HB-1100**    **OA-DA-3006**



An OxyArm™ Product

*1st Open Oxygen System™*  
*We are the experts!*

Patented

OXYMASK IS A REGISTERED TRADEMARK



705-726-9383  
 email: [smedic@southmedic.com](mailto:smedic@southmedic.com)

[www.southmedic.com](http://www.southmedic.com)    [www.oxyarm.com](http://www.oxyarm.com)  
 50 Alliance Blvd., Barrie, ON., Canada L4M 5K3  
 ISO 13485