

OxyMask ETCO₂

**DELIVERS 24% to 65% OXYGEN
with CO₂ Monitoring**



OxyMask™

OXYMASK

DELIVERS 24% to 65% OXYGEN with Endtidal CO₂

OxyMask's & OxyChin's logical advantages:

- Eliminates CO₂ rebreathing
- Provides quick access for oral suctioning
- Ideal for conscious sedation & surgical post-op oxygenation
- Open design allows communication between patient & caregiver
- OxyMask delivers 1-15 LPM (with CO₂ Monitoring)
- Reduces claustrophobia
- Safe when emesis occurs

Additional OxyChin advantages:

- Ideal for eye, ear, nose or neuro surgery
- OxyChin delivers 1-15 LPM (with CO₂ monitoring)

Safe Effective Monitoring of CO₂



OxyMask
Adult

OxyPlus
25% Larger

Over Ear
Dual Elastics

OxyChin

OxyKid
Endtidal

ORDER CODES:

- OM-2125-8 OxyMask Endtidal CO₂ - box of 25
- OM-2425-8 OxyMask Endtidal CO₂ with dual ear elastics - box of 25
- OP-2125-8 OxyPlus Endtidal CO₂ - 25% larger than OxyMask - box of 25
- OP-2425-8 OxyPlus Endtidal CO₂ with dual ear elastics - box of 25
- OC-2025 OxyChin Endtidal CO₂ without tubing - box of 25
- OC-2025-8 OxyChin Endtidal CO₂ with tubing - box of 25
- OK-2125-8 OxyKid Endtidal CO₂ with tubing - box of 25

OxyMask[™]
An OxyArm[™] Product

*1st Open Oxygen System.
We are the experts!*
Patented

Oxygen Flow Settings	Approximate O ₂ Concentrations
1 LPM	24% - 25%
2 LPM	27% - 32%
4 LPM	34% - 40%
6 LPM	42% - 48%
8 LPM	49% - 55%
10 LPM	53% - 58%
12-15 LPM	59% - 65%



CO₂ Luer
Connection

Southmedic

705-726-9383
email: smedic@southmedic.com

www.southmedic.com www.oxyarm.com

50 Alliance Blvd., Barrie, ON., Canada L4M 5K3
ISO 13485