

Sterilization Tray

NEW



- colour coded for procedural differentiation
- compatible with steam or Sterrad™ sterilization
- easy to open and close
- interchangeable tops and bottoms
- clear and strong construction
- used for tracking, storage and sterilization
- one year or 500 cycle warranty

Sterilization Tray

- Colour coding allows for easy department differentiation
- Larger holes for quicker liquid evaporation and steam penetration
- Base construction keeps silicone pin mat raised, facilitating air flow and eliminating pooling of liquid
- Easy open / close latch eliminates torn gloves
- Lids and bases are engineered identically, eliminating mis-matched components
- Durable construction - one year or 500 cycle warranty
- Transparent design for easy viewing of contents
- Multi-purpose containers - storage/ protection
 - tracking
 - sterilization

Ideal for:

- Small boxes
 - assorted small parts
 - burrs / drill bits
 - cardio & ophthalmic instruments
- Medium + large boxes
 - micro instruments
 - ophthalmic, neuro & plastics
 - dental instruments



Burr block available

Products Codes:

	TRAY WITH MAT	INNER DIMENSIONS
STS-RE	Sterilization Tray, small - Red	1.75"W x 2.75" L x 1"
STS-GR	Sterilization Tray, small - Green	1.75"W x 2.75" L x 1"
STS-NA	Sterilization Tray, small - Natural	1.75"W x 2.75" L x 1"
STM-RE	Sterilization Tray, medium - Red	3"W x 6" L x 1"
STM-GR	Sterilization Tray, medium - Green	3"W x 6" L x 1"
STM-NA	Sterilization Tray, medium - Natural	3"W x 6" L x 1"
STL-RE	Sterilization Tray, large - Red	3"W x 7 3/8" L x 1"
STL-GR	Sterilization Tray, large - Green	3"W x 7 3/8" L x 1"
STL-NA	Sterilization Tray, large - Natural	3"W x 7 3/8" L x 1"

Box colours available:



	TRAY WITHOUT MAT	
STS-RE-NM	Sterilization Tray, small - Red	1.75"W x 2.75" L x 1"
STS-GR-NM	Sterilization Tray, small - Green	1.75"W x 2.75" L x 1"
STS-NA-NM	Sterilization Tray, small - Natural	1.75"W x 2.75" L x 1"
STM-RE-NM	Sterilization Tray, medium - Red	3"W x 6" L x 1"
STM-GR-NM	Sterilization Tray, medium - Green	3"W x 6" L x 1"
STM-NA-NM	Sterilization Tray, medium - Natural	3"W x 6" L x 1"
STL-RE-NM	Sterilization Tray, large - Red	3"W x 7 3/8" L x 1"
STL-GR-NM	Sterilization Tray, large - Green	3"W x 7 3/8" L x 1"
STL-NA-NM	Sterilization Tray, large - Natural	3"W x 7 3/8" L x 1"

SBB-1000	Sterilization Burr Block	
MAT1000	Pin Cushion Mat	13"x18"

(Large pin cushion mats available - allowing you to cut to your specific requirements)



50 Alliance Blvd., Barrie, ON., Canada L4M 5K3
 705-726-9383 • 705-726-9537 Fax ISO 13485
 Email: smedic@southmedic.on.ca

EDI 406000117

www.oxyarm.com • www.southmedic.com