

# Introducing The New Vapofil



## With Stronger and More Reliable Features

### Benefits

- Faster filling
- More flexible / No kinking
- Improved flow
- No air locks
- Improved seal / Less waste
- Absolute dimensional stability
- Polished stainless steel prevents scratching of the vaporizer gasket
- Enhanced sealing at block end
  - Eliminate leakage
- Enhanced product safety

# Vaporizer Filler Tubes That Work First Drop to Last

## Vapofil™



### Steel Ordering Info

Vaporizer	Cat. No.
■ Florane/Isoflurane	V0507F
■ Ultane/Sevoflurane	V0507S

### Nylon Ordering Info

Vaporizer	Cat. No.
■ Florane/Isoflurane	8907-F
■ Halothane/Fluothane	8907-H
■ Ultane/Sevoflurane	8907-S

## AntiSpil™



### Nylon Ordering Info

Vaporizer	Cat. No.
■ Florane/Isoflurane	9010-F
■ Halothane/Fluothane	9010-H
■ Ultane/Sevoflurane	9010-S

## Features of Steel:

- Convoluted outer tubing
- Inner tube corrugated at bend point
- Colour coded and keyed for agent specific use
- Separate air and fluid tube
- Anti-pollution plug
- Rapid filling
- Durable
- Holder available



## Features of Nylon:

- Flexible outer tubing
- Separate inner tubes for air and agent
- Inner sleeve reinforced by brass tubing
- Colour coded and keyed for agent specific use
- Anti-pollution plug
- Rapid filling
- Durable
- Holder available



## Features:

- Curved tip and flow meter let you pour the agent with less risk of spills
- Colour coded and keyed cap for agent specific use
- Anti-pollution plug
- Durable
- Holder available



50 Alliance Blvd., Barrie, Ontario, Canada L4M 5K3

1-705-726-9383

1-705-728-9537

FAX [www.southmedic.com](http://www.southmedic.com)

Email: [smedic@southmedic.on.ca](mailto:smedic@southmedic.on.ca)

EDI 406000117

CE marked (Europe) • Issued substantial equivalence from FDA (USA) • Licensed under Health Canada

11.06 Ver1