

NeVAP Secure-Swap Endotracheal Tube Securement System

Smarter Securement for Safer, Cost-Efficient Airway Care

Distributed Exclusively in Canada by Southmedic

Overview

Designed for **ICU and other critical-care environments**, the NeVAP Secure-Swap ETT device provides advanced securement for **orally intubated patients**, replacing traditional taping methods and other problematic securement devices, with a comfortable, adjustable interface that protects skin integrity and simplifies tube repositioning.

Its air-inflated cheek cuffs disperse pressure evenly, allowing clinicians to monitor facial skin and prevent device-related pressure injuries—without interrupting ventilation or stability.

The Secure-Swap Endotracheal Tube Holder allows clinicians to practice **skin breakdown prevention protocols that align with National Pressure Ulcer Advisory Panel (NPUAP) recommendations and guidelines**, helping reduce the risk of medical device-related pressure injuries (MDRPI) and improve overall patient outcomes.

When components become soiled or require change, only the attachment is replaced—eliminating the need to discard the entire device.

Clinical Applications

- For orally intubated or ventilated patients in the **ICU, OR, and ED units**.
- Ideal for patients requiring prolonged intubation or frequent repositioning.
- Not typically used in short-term emergency intubations.

Key Features & Benefits

Protects Skin Integrity

- Air-inflated cheek cuffs distribute pressure to reduce skin breakdown.
- Clear, open-center cuffs allow continuous visibility of the patient's skin for early detection of redness or irritation.
- Heat and moisture are vented away from the face, promoting comfort and compliance.

Improves Clinical Efficiency & Safety

- Secure, adjustable design replaces face taping for a consistent seal and better airway management.
- Allows side-to-side repositioning to relieve pressure points during long-term ventilation.
- Gentle adhesive backing on the cheek cuffs ensures safe securement.
- Compatible with a wide range of endotracheal tubes and configurations—including options with or without bite block.
- Accessory pieces detach easily, so the full system never needs to be replaced mid-care.

Reduces Cost & Waste

- Replace only what's required—each component can be swapped without replacing the frame.
- Modular design supports infection-control policies and extends device life between changes.



Modular Attachments for Smarter, Cost-Effective Care

The Secure-Swap device features a modular design that lets clinicians replace only what's needed—saving time, cost, and waste. Each attachment connects securely to the base frame (Secure-Swap F), supporting various airway management needs without replacing the entire device.

- **ETF102 — Secure-Swap F (Base Frame)** The foundation of the system. Attachments connect directly to this reusable base, allowing fast swaps when components become soiled or protocols require change.
- **ETF101 — Secure-Swap G** For non-Medtronic/Covidien tubes. Includes a built-in bite block and anti-kink design to prevent unplanned extubation while supporting the natural McGill curve of the ETT.
- **ETF103 — Secure-Swap O** Also for non-Medtronic/Covidien tubes. Offers the same anti-kink and migration-prevention features as the G but without a bite block for flexibility in patient care.
- **ETF104 — Secure-Swap A** A universal securement option that stabilizes the body of the ETT and other devices such as orogastric tubes and pressure monitoring probes.
- **ETF105 — Secure-Swap B** A universal attachment with a straight bite block, ideal when added protection against occlusion or biting is needed.



ETF1022 — Secure-Swap F (ETF1052 — Secure-Swap B attached)



ETF1012 — Secure-Swap G



ETF1032 — Secure-Swap O



ETF1042 — Secure-Swap A



ETF1052 — Secure-Swap B

PART NUMBER	PRODUCT DESCRIPTION	QUANTITY
ETF-102	Secure-Swap F - Base Frame	20/box
ETF-104	Secure-Swap A - ETT holder attachment	20/box
ETF-105	Secure-Swap B - ETT holder attachment with bite guard	20/box
ETF-101	Secure-Swap G - ETT holder sleeve with bite guard	20/box
ETF-103	Secure-Swap O - ETT holder sleeve	20/box

Finding a Better Way in Airway Care

At Southmedic, we're proud to bring innovative solutions like NeVAP Secure-Swap to healthcare professionals across Canada—helping improve patient safety, protect skin integrity, and reduce costs through smarter securement.

Talk to your Southmedic representative to learn more or request a product demonstration.



Distributed by Southmedic
50 Alliance Blvd.,
Barrie, ON Canada L4M 5K3
www.southmedic.com

Email: custserv@southmedic.com
1-800-463-7146 Toll free in North America
1-705-720-1902
1-705-728-9537 Fax

