


# Medical Grade Silicone Tubing

Not Made with Natural Rubber Latex



## Engineered for Purity and Safety

Platinum-cured for the highest degree of purity, this medical grade silicone tubing meets USP Class VI requirements and is not made with BPA, phthalates or natural rubber latex. It is naturally translucent allowing for easy flow visibility and its elasticity supports easy movement.

## Why Choose Medical Grade Silicone Tubing

- **High purity and safe** – platinum cured and not made with BPA, phthalates or natural rubber latex
- **Biocompatible and compliant** – meets USP Class VI; EU RoHS 3 and EU REACH SVHC compliant
- **Clear visibility** – naturally translucent for quick flow inspection and monitoring
- **Durable and flexible** – highly elastic and made without phthalates.
- **Sterilization-ready** – supplied in individually heat-sealed polybags

## Cleaning and Sterilization Recommendations

To obtain cleaning and sterilization recommendations, please reach out to your local territory representative.

### Ordering Information

PRODUCT CODE	DESCRIPTION
T2015-100	Silicone Tubing .188" x .313" x 100'
T2018-100	Silicone Tubing .250" x .438" x 100'
T2016-100	Silicone Tubing .188" x .375" x 100'
T2017-100	Silicone Tubing .250" x .375" x 100'
T2013-100	Silicone Tubing .125" x .250" x 100'
T2021-100	Silicone Tubing .313" x .500" x 100'
T2019-100	Silicone Tubing .250" x .500" x 100'
T2020-100	Silicone Tubing .313" x .438" x 100'
T2022-100	Silicone Tubing .375" x .500" x 100'

