


# DynaStretch®

Clinical and cosmetic outcomes matter


## Pre-surgical Skin Expansion

DynaStretch strips are designed to aid in pre-surgical skin expansion. This expansion occurs secondarily to the continued cyclic stretching of the skin by Dynamic Tissue Systems. New surplus skin may be consistent with the surrounding tissue and can be used to span a tissue defect created by surgical excision. Gentle, dynamic tissue stretching prior to planned excision avoids leaving a skin defect and allows for a sound primary closure.

This system is non-invasive and can be applied days or weeks pre-op without hindering patient mobility or quality of life. Uses may include any procedure in which a skin graft site expansion is desired.



### Results at a Glance\*



DynaStretch used to pre-operatively stretch the area surrounding the skin lesion

DSX24 DynaStretch (0.5 x 4.5")

Small DynaStretch Strips (4/pack) Qty: 60 (15 packs/set)

DSX38 DynaStretch (0.75 x 8.5")

Large DynaStretch Strips (4/pack) Qty: 60 (15 packs/set)

\* This case contains the opinions of and personal techniques practiced by the treating physician. The techniques presented herein are for informational purposes only. The decision of which techniques to use in a particular clinical application lies with the physician based on patient profile, particular circumstances surrounding the procedure, and previous clinical experiences.

**Southmedic**

Dynamic Tissue Systems

Respiratory

Surgical

Blades & Scalpels

Southmedic Inc.,  
50 Alliance Blvd.,  
Barrie, ON  
Canada L4M 5K3

www.dynamicissuesystems.com  
705-726-9383  
1-800-463-7146 Toll free in North America  
contactus@southmedic.com